

Table of contents				
\\USER				
AP				
TEST				
EMPRISE_220217				
QuinPilot_32Ch				
cmrr_mbep2d_bold_2.2iso_mod270122				
cmrr_mbep2d_bold_2.2iso_mod270122				
cmrr_mbep2d_bold_2.2iso_mod270122				
cmrr_mbep2d_bold_2.2iso_mod270122				
cmrr_mbep2d_bold_2.2iso_mod270122				
cmrr_mbep2d_bold_2.2iso_mod270122				
cmrr_mbep2d_bold_2.2iso_mod270122_loc_word				
cmrr_mbep2d_bold_2.2iso_mod270122_loc_visual				

\\USER\AP\TEST\EMPRISE\_220217\QuinPilot\_32Ch

TA: 0:47 PM: REF Voxel size: 1.0×1.0×5.0 mmPAT: Off Rel. SNR: 1.00 : tfl

**Properties**

Prio recon	Off
Load images to viewer	On
Inline movie	Off
Auto store images	On
Load images to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Auto close inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	Single measurement

**Routine**

Slice group	1
Slices	3
Dist. factor	50 %
Position	L0.0 P13.5 H0.0 mm
Orientation	Sagittal
Phase enc. dir.	A >> P
Slice group	2
Slices	1
Dist. factor	50 %
Position	L0.0 P2.8 H0.0 mm
Orientation	Transversal
Phase enc. dir.	A >> P
Slice group	3
Slices	1
Dist. factor	50 %
Position	L0.0 P12.6 H3.7 mm
Orientation	Coronal
Phase enc. dir.	R >> L
AutoAlign	---
Phase oversampling	0 %
FoV read	250 mm
FoV phase	100.0 %
Slice thickness	5.0 mm
TR	3000.0 ms
TE	3.63 ms
Averages	3
Concatenations	5
Filter	Prescan Normalize
Coil elements	HEA;HEP

**Contrast - Common**

TR	3000.0 ms
TE	3.63 ms
TD	0 ms
Magn. preparation	Slice-sel. IR
TI	1500 ms
Flip angle	10 deg
Fat suppr.	None
Water suppr.	None

**Contrast - Dynamic**

Averages	3
Averaging mode	Short term
Reconstruction	Magnitude
Measurements	1
Multiple series	Each measurement

**Resolution - Common**

FoV read	250 mm
FoV phase	100.0 %
Slice thickness	5.0 mm
Base resolution	256
Phase resolution	100 %
Phase partial Fourier	Off
Interpolation	Off

**Resolution - iPAT**

PAT mode	None
----------	------

**Resolution - Filter Image**

Image Filter	Off
Distortion Corr.	Off
Prescan Normalize	On
Unfiltered images	Off
Normalize	Off
B1 filter	Off

**Resolution - Filter Rawdata**

Raw filter	Off
Elliptical filter	Off

**Geometry - Common**

Slice group	1
Slices	3
Dist. factor	50 %
Position	L0.0 P13.5 H0.0 mm
Orientation	Sagittal
Phase enc. dir.	A >> P
Slice group	2
Slices	1
Dist. factor	50 %
Position	L0.0 P2.8 H0.0 mm
Orientation	Transversal
Phase enc. dir.	A >> P
Slice group	3
Slices	1
Dist. factor	50 %
Position	L0.0 P12.6 H3.7 mm
Orientation	Coronal
Phase enc. dir.	R >> L
FoV read	250 mm
FoV phase	100.0 %
Slice thickness	5.0 mm
TR	3000.0 ms
Multi-slice mode	Sequential
Series	Interleaved
Concatenations	5

**Geometry - AutoAlign**

Slice group	1
Position	L0.0 P13.5 H0.0 mm
Orientation	Sagittal
Phase enc. dir.	A >> P
Slice group	2
Position	L0.0 P2.8 H0.0 mm
Orientation	Transversal
Phase enc. dir.	A >> P
Slice group	3

**Geometry - AutoAlign**

Position	L0.0 P12.6 H3.7 mm
Orientation	Coronal
Phase enc. dir.	R >> L
AutoAlign	---
Initial Position	L0.0 P13.5 H0.0
L	0.0 mm
P	13.5 mm
H	0.0 mm
Initial Rotation	0.00 deg
Initial Orientation	Sagittal

**Geometry - Navigator****System - Miscellaneous**

Positioning mode	REF
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Adaptive Combine
Save uncombined	Off
Matrix Optimization	Off
AutoAlign	---
Coil Select Mode	On - AutoCoilSelect

**System - Adjustments**

B0 Shim mode	Tune up
B1 Shim mode	TrueForm
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto

**System - Adjust Volume**

Position	Isocenter
Orientation	Transversal
Rotation	0.00 deg
A >> P	263 mm
R >> L	350 mm
F >> H	350 mm
Reset	Off

**System - Tx/Rx**

Frequency 1H	123.253584 MHz
Correction factor	1
Gain	High
Img. Scale Cor.	1.000
Reset	Off
? Ref. amplitude 1H	0.000 V

**Physio - Signal1**

1st Signal/Mode	None
TR	3000.0 ms
Concatenations	5

**Physio - Cardiac**

Magn. preparation	Slice-sel. IR
TI	1500 ms
Fat suppr.	None
Dark blood	Off

**Physio - Cardiac**

FoV read	250 mm
FoV phase	100.0 %
Phase resolution	100 %

**Physio - PACE**

Resp. control	Off
Concatenations	5

**Inline - Common**

Subtract	Off
Measurements	1
StdDev	Off
Save original images	On

**Inline - MIP**

MIP-Sag	Off
MIP-Cor	Off
MIP-Tra	Off
MIP-Time	Off
Save original images	On

**Inline - Composing**

Distortion Corr.	Off
------------------	-----

**Inline - Maplt**

Save original images	On
Maplt	None
Flip angle	10 deg
Measurements	1
TR	3000.0 ms
TE	3.63 ms

**Sequence - Part 1**

Introduction	On
Dimension	2D
Asymmetric echo	Off
Flow comp.	No
Multi-slice mode	Sequential
Echo spacing	7.1 ms
Bandwidth	240 Hz/Px

**Sequence - Part 2**

RF pulse type	Normal
Gradient mode	Whisper
Excitation	Slice-sel.
RF spoiling	On
Incr. Gradient spoiling	Off
Turbo factor	256

**Sequence - Assistant**

Mode	Off
------	-----

\\USER\AP\TEST\EMPRISE\_220217\cmrr\_mbep2d\_bold\_2.2iso\_mod270122

TA: 6:06 PM: FIX Voxel size: 2.2×2.2×2.2 mmPAT: Off Rel. SNR: 1.00 : efpid

**Properties**

Prio recon	Off
Load images to viewer	On
Inline movie	Off
Auto store images	On
Load images to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Auto close inline display	Off
Start measurement without further preparation	Off
Wait for user to start	On
Start measurements	Single measurement

**Routine**

Slice group	1
Slices	51
Dist. factor	13 %
Position	R1.9 P25.7 H32.9 mm
Orientation	T > C-15.4
Phase enc. dir.	A >> P
AutoAlign	---
Phase oversampling	0 %
FoV read	204 mm
FoV phase	100.0 %
Slice thickness	2.20 mm
TR	2000 ms
TE	24.60 ms
Multi-band accel. factor	3
Filter	Prescan Normalize
Coil elements	HEA;HEP

**Contrast - Common**

TR	2000 ms
TE	24.60 ms
MTC	Off
Magn. preparation	None
Flip angle	80 deg
Fat suppr.	Fat sat.

**Contrast - Dynamic**

Averaging mode	Long term
Reconstruction	Magnitude
Measurements	178
Delay in TR	0 ms
Multiple series	Off

**Resolution - Common**

FoV read	204 mm
FoV phase	100.0 %
Slice thickness	2.20 mm
Base resolution	92
Phase resolution	100 %
Phase partial Fourier	6/8
Interpolation	Off

**Resolution - iPAT**

PAT mode	None
----------	------

**Resolution - Filter Image**

Distortion Corr.	Off
------------------	-----

**Resolution - Filter Image**

Prescan Normalize	On
-------------------	----

**Resolution - Filter Rawdata**

Raw filter	Off
Elliptical filter	Off
Hamming	Off

**Geometry - Common**

Slice group	1
Slices	51
Dist. factor	13 %
Position	R1.9 P25.7 H32.9 mm
Orientation	T > C-15.4
Phase enc. dir.	A >> P
FoV read	204 mm
FoV phase	100.0 %
Slice thickness	2.20 mm
TR	2000 ms
Multi-slice mode	Interleaved
Series	Interleaved
Multi-band accel. factor	3

**Geometry - AutoAlign**

Slice group	1
Position	R1.9 P25.7 H32.9 mm
Orientation	T > C-15.4
Phase enc. dir.	A >> P
AutoAlign	---
Initial Position	R1.9 P25.7 H32.9
R	1.9 mm
P	25.7 mm
H	32.9 mm
Initial Rotation	0.00 deg
Initial Orientation	T > C
T > C	-15.4
> S	0.0

**Geometry - Saturation**

Fat suppr.	Fat sat.
Special sat.	None

**System - Miscellaneous**

Positioning mode	FIX
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Sum of Squares
Matrix Optimization	Off
AutoAlign	---
Coil Select Mode	Default

**System - Adjustments**

B0 Shim mode	Standard
B1 Shim mode	TrueForm
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Dominant Fat	Off

**System - Adjustments**

Assume Silicone	Off
Adjustment Tolerance	Auto

**System - Adjust Volume**

Position	R1.9 P25.7 H32.9 mm
Orientation	T > C-15.4
Rotation	0.00 deg
A >> P	204 mm
R >> L	204 mm
F >> H	127 mm
Reset	Off

**System - Tx/Rx**

Frequency 1H	123.253584 MHz
Correction factor	1
Gain	High
Img. Scale Cor.	1.000
Reset	Off
? Ref. amplitude 1H	0.000 V

**Physio - Signal1**

1st Signal/Mode	None
TR	2000 ms
Multi-band accel. factor	3

**BOLD**

GLM Statistics	Off
Dynamic t-maps	Off
Ignore meas. at start	0
Ignore after transition	0
Model transition states	On
Temp. highpass filter	On
Threshold	4.00
Paradigm size	20
Meas[1]	Baseline
Meas[2]	Baseline
Meas[3]	Baseline
Meas[4]	Baseline
Meas[5]	Baseline
Meas[6]	Baseline
Meas[7]	Baseline
Meas[8]	Baseline
Meas[9]	Baseline
Meas[10]	Baseline
Meas[11]	Active
Meas[12]	Active
Meas[13]	Active
Meas[14]	Active
Meas[15]	Active
Meas[16]	Active
Meas[17]	Active
Meas[18]	Active
Meas[19]	Active
Meas[20]	Active
Motion correction	Off
Spatial filter	Off
Measurements	178
Delay in TR	0 ms
Multiple series	Off

**Sequence - Part 1**

Introduction	Off
Contrasts	1
Flow comp.	No

**Sequence - Part 1**

Multi-slice mode	Interleaved
Free echo spacing	Off
Echo spacing	0.77 ms
Bandwidth	1552 Hz/Px

**Sequence - Part 2**

EPI factor	92
Gradient mode	Normal
RF spoiling	Off

**Sequence - Special**

Excite pulse duration	4380 us
Single-band images	Off
MB LeakBlock kernel	On
MB dual kernel	Off
MB RF phase scramble	Off
SENSE1 coil combine	On
Invert RO/PE polarity	Off
PF omits higher k-space	Off
Disable freq. update	Off
Force equal slice timing	Off
Online multi-band recon.	Online
FFT scale factor	1.00
Physio recording	Off
Triggering scheme	Standard

\\USER\AP\TEST\EMPRISE\_220217\cmrr\_mbep2d\_bold\_2.2iso\_mod270122

TA: 6:06 PM: FIX Voxel size: 2.2×2.2×2.2 mmPAT: Off Rel. SNR: 1.00 : efpid

**Properties**

Prio recon	Off
Load images to viewer	On
Inline movie	Off
Auto store images	On
Load images to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Auto close inline display	Off
Start measurement without further preparation	Off
Wait for user to start	On
Start measurements	Single measurement

**Routine**

Slice group	1
Slices	51
Dist. factor	13 %
Position	R1.9 P25.7 H32.9 mm
Orientation	T > C-15.4
Phase enc. dir.	A >> P
AutoAlign	---
Phase oversampling	0 %
FoV read	204 mm
FoV phase	100.0 %
Slice thickness	2.20 mm
TR	2000 ms
TE	24.60 ms
Multi-band accel. factor	3
Filter	Prescan Normalize
Coil elements	HEA;HEP

**Contrast - Common**

TR	2000 ms
TE	24.60 ms
MTC	Off
Magn. preparation	None
Flip angle	80 deg
Fat suppr.	Fat sat.

**Contrast - Dynamic**

Averaging mode	Long term
Reconstruction	Magnitude
Measurements	178
Delay in TR	0 ms
Multiple series	Off

**Resolution - Common**

FoV read	204 mm
FoV phase	100.0 %
Slice thickness	2.20 mm
Base resolution	92
Phase resolution	100 %
Phase partial Fourier	6/8
Interpolation	Off

**Resolution - iPAT**

PAT mode	None
----------	------

**Resolution - Filter Image**

Distortion Corr.	Off
------------------	-----

**Resolution - Filter Image**

Prescan Normalize	On
-------------------	----

**Resolution - Filter Rawdata**

Raw filter	Off
Elliptical filter	Off
Hamming	Off

**Geometry - Common**

Slice group	1
Slices	51
Dist. factor	13 %
Position	R1.9 P25.7 H32.9 mm
Orientation	T > C-15.4
Phase enc. dir.	A >> P
FoV read	204 mm
FoV phase	100.0 %
Slice thickness	2.20 mm
TR	2000 ms
Multi-slice mode	Interleaved
Series	Interleaved
Multi-band accel. factor	3

**Geometry - AutoAlign**

Slice group	1
Position	R1.9 P25.7 H32.9 mm
Orientation	T > C-15.4
Phase enc. dir.	A >> P
AutoAlign	---
Initial Position	R1.9 P25.7 H32.9
R	1.9 mm
P	25.7 mm
H	32.9 mm
Initial Rotation	0.00 deg
Initial Orientation	T > C
T > C	-15.4
> S	0.0

**Geometry - Saturation**

Fat suppr.	Fat sat.
Special sat.	None

**System - Miscellaneous**

Positioning mode	FIX
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Sum of Squares
Matrix Optimization	Off
AutoAlign	---
Coil Select Mode	Default

**System - Adjustments**

B0 Shim mode	Standard
B1 Shim mode	TrueForm
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Dominant Fat	Off

**System - Adjustments**

Assume Silicone	Off
Adjustment Tolerance	Auto

**System - Adjust Volume**

Position	R1.9 P25.7 H32.9 mm
Orientation	T > C-15.4
Rotation	0.00 deg
A >> P	204 mm
R >> L	204 mm
F >> H	127 mm
Reset	Off

**System - Tx/Rx**

Frequency 1H	123.253584 MHz
Correction factor	1
Gain	High
Img. Scale Cor.	1.000
Reset	Off
? Ref. amplitude 1H	0.000 V

**Physio - Signal1**

1st Signal/Mode	None
TR	2000 ms
Multi-band accel. factor	3

**BOLD**

GLM Statistics	Off
Dynamic t-maps	Off
Ignore meas. at start	0
Ignore after transition	0
Model transition states	On
Temp. highpass filter	On
Threshold	4.00
Paradigm size	20
Meas[1]	Baseline
Meas[2]	Baseline
Meas[3]	Baseline
Meas[4]	Baseline
Meas[5]	Baseline
Meas[6]	Baseline
Meas[7]	Baseline
Meas[8]	Baseline
Meas[9]	Baseline
Meas[10]	Baseline
Meas[11]	Active
Meas[12]	Active
Meas[13]	Active
Meas[14]	Active
Meas[15]	Active
Meas[16]	Active
Meas[17]	Active
Meas[18]	Active
Meas[19]	Active
Meas[20]	Active
Motion correction	Off
Spatial filter	Off
Measurements	178
Delay in TR	0 ms
Multiple series	Off

**Sequence - Part 1**

Introduction	Off
Contrasts	1
Flow comp.	No

**Sequence - Part 1**

Multi-slice mode	Interleaved
Free echo spacing	Off
Echo spacing	0.77 ms
Bandwidth	1552 Hz/Px

**Sequence - Part 2**

EPI factor	92
Gradient mode	Normal
RF spoiling	Off

**Sequence - Special**

Excite pulse duration	4380 us
Single-band images	Off
MB LeakBlock kernel	On
MB dual kernel	Off
MB RF phase scramble	Off
SENSE1 coil combine	On
Invert RO/PE polarity	Off
PF omits higher k-space	Off
Disable freq. update	Off
Force equal slice timing	Off
Online multi-band recon.	Online
FFT scale factor	1.00
Physio recording	Off
Triggering scheme	Standard

\\USER\AP\TEST\EMPRISE\_220217\cmrr\_mbep2d\_bold\_2.2iso\_mod270122

TA: 6:06 PM: FIX Voxel size: 2.2×2.2×2.2 mmPAT: Off Rel. SNR: 1.00 : efpid

**Properties**

Prio recon	Off
Load images to viewer	On
Inline movie	Off
Auto store images	On
Load images to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Auto close inline display	Off
Start measurement without further preparation	Off
Wait for user to start	On
Start measurements	Single measurement

**Routine**

Slice group	1
Slices	51
Dist. factor	13 %
Position	R1.9 P25.7 H32.9 mm
Orientation	T > C-15.4
Phase enc. dir.	A >> P
AutoAlign	---
Phase oversampling	0 %
FoV read	204 mm
FoV phase	100.0 %
Slice thickness	2.20 mm
TR	2000 ms
TE	24.60 ms
Multi-band accel. factor	3
Filter	Prescan Normalize
Coil elements	HEA;HEP

**Contrast - Common**

TR	2000 ms
TE	24.60 ms
MTC	Off
Magn. preparation	None
Flip angle	80 deg
Fat suppr.	Fat sat.

**Contrast - Dynamic**

Averaging mode	Long term
Reconstruction	Magnitude
Measurements	178
Delay in TR	0 ms
Multiple series	Off

**Resolution - Common**

FoV read	204 mm
FoV phase	100.0 %
Slice thickness	2.20 mm
Base resolution	92
Phase resolution	100 %
Phase partial Fourier	6/8
Interpolation	Off

**Resolution - iPAT**

PAT mode	None
----------	------

**Resolution - Filter Image**

Distortion Corr.	Off
------------------	-----

**Resolution - Filter Image**

Prescan Normalize	On
-------------------	----

**Resolution - Filter Rawdata**

Raw filter	Off
Elliptical filter	Off
Hamming	Off

**Geometry - Common**

Slice group	1
Slices	51
Dist. factor	13 %
Position	R1.9 P25.7 H32.9 mm
Orientation	T > C-15.4
Phase enc. dir.	A >> P
FoV read	204 mm
FoV phase	100.0 %
Slice thickness	2.20 mm
TR	2000 ms
Multi-slice mode	Interleaved
Series	Interleaved
Multi-band accel. factor	3

**Geometry - AutoAlign**

Slice group	1
Position	R1.9 P25.7 H32.9 mm
Orientation	T > C-15.4
Phase enc. dir.	A >> P
AutoAlign	---
Initial Position	R1.9 P25.7 H32.9
R	1.9 mm
P	25.7 mm
H	32.9 mm
Initial Rotation	0.00 deg
Initial Orientation	T > C
T > C	-15.4
> S	0.0

**Geometry - Saturation**

Fat suppr.	Fat sat.
Special sat.	None

**System - Miscellaneous**

Positioning mode	FIX
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Sum of Squares
Matrix Optimization	Off
AutoAlign	---
Coil Select Mode	Default

**System - Adjustments**

B0 Shim mode	Standard
B1 Shim mode	TrueForm
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Dominant Fat	Off



**System - Adjustments**

Assume Silicone	Off
Adjustment Tolerance	Auto

**System - Adjust Volume**

Position	R1.9 P25.7 H32.9 mm
Orientation	T > C-15.4
Rotation	0.00 deg
A >> P	204 mm
R >> L	204 mm
F >> H	127 mm
Reset	Off

**System - Tx/Rx**

Frequency 1H	123.253584 MHz
Correction factor	1
Gain	High
Img. Scale Cor.	1.000
Reset	Off
? Ref. amplitude 1H	0.000 V

**Physio - Signal1**

1st Signal/Mode	None
TR	2000 ms
Multi-band accel. factor	3

**BOLD**

GLM Statistics	Off
Dynamic t-maps	Off
Ignore meas. at start	0
Ignore after transition	0
Model transition states	On
Temp. highpass filter	On
Threshold	4.00
Paradigm size	20
Meas[1]	Baseline
Meas[2]	Baseline
Meas[3]	Baseline
Meas[4]	Baseline
Meas[5]	Baseline
Meas[6]	Baseline
Meas[7]	Baseline
Meas[8]	Baseline
Meas[9]	Baseline
Meas[10]	Baseline
Meas[11]	Active
Meas[12]	Active
Meas[13]	Active
Meas[14]	Active
Meas[15]	Active
Meas[16]	Active
Meas[17]	Active
Meas[18]	Active
Meas[19]	Active
Meas[20]	Active
Motion correction	Off
Spatial filter	Off
Measurements	178
Delay in TR	0 ms
Multiple series	Off

**Sequence - Part 1**

Introduction	Off
Contrasts	1
Flow comp.	No

**Sequence - Part 1**

Multi-slice mode	Interleaved
Free echo spacing	Off
Echo spacing	0.77 ms
Bandwidth	1552 Hz/Px

**Sequence - Part 2**

EPI factor	92
Gradient mode	Normal
RF spoiling	Off

**Sequence - Special**

Excite pulse duration	4380 us
Single-band images	Off
MB LeakBlock kernel	On
MB dual kernel	Off
MB RF phase scramble	Off
SENSE1 coil combine	On
Invert RO/PE polarity	Off
PF omits higher k-space	Off
Disable freq. update	Off
Force equal slice timing	Off
Online multi-band recon.	Online
FFT scale factor	1.00
Physio recording	Off
Triggering scheme	Standard

\\USER\AP\TEST\EMPRISE\_220217\cmrr\_mbep2d\_bold\_2.2iso\_mod270122

TA: 6:06 PM: FIX Voxel size: 2.2×2.2×2.2 mmPAT: Off Rel. SNR: 1.00 : efpid

**Properties**

Prio recon	Off
Load images to viewer	On
Inline movie	Off
Auto store images	On
Load images to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Auto close inline display	Off
Start measurement without further preparation	Off
Wait for user to start	On
Start measurements	Single measurement

**Routine**

Slice group	1
Slices	51
Dist. factor	13 %
Position	R1.9 P25.7 H32.9 mm
Orientation	T > C-15.4
Phase enc. dir.	A >> P
AutoAlign	---
Phase oversampling	0 %
FoV read	204 mm
FoV phase	100.0 %
Slice thickness	2.20 mm
TR	2000 ms
TE	24.60 ms
Multi-band accel. factor	3
Filter	Prescan Normalize
Coil elements	HEA;HEP

**Contrast - Common**

TR	2000 ms
TE	24.60 ms
MTC	Off
Magn. preparation	None
Flip angle	80 deg
Fat suppr.	Fat sat.

**Contrast - Dynamic**

Averaging mode	Long term
Reconstruction	Magnitude
Measurements	178
Delay in TR	0 ms
Multiple series	Off

**Resolution - Common**

FoV read	204 mm
FoV phase	100.0 %
Slice thickness	2.20 mm
Base resolution	92
Phase resolution	100 %
Phase partial Fourier	6/8
Interpolation	Off

**Resolution - iPAT**

PAT mode	None
----------	------

**Resolution - Filter Image**

Distortion Corr.	Off
------------------	-----

**Resolution - Filter Image**

Prescan Normalize	On
-------------------	----

**Resolution - Filter Rawdata**

Raw filter	Off
Elliptical filter	Off
Hamming	Off

**Geometry - Common**

Slice group	1
Slices	51
Dist. factor	13 %
Position	R1.9 P25.7 H32.9 mm
Orientation	T > C-15.4
Phase enc. dir.	A >> P
FoV read	204 mm
FoV phase	100.0 %
Slice thickness	2.20 mm
TR	2000 ms
Multi-slice mode	Interleaved
Series	Interleaved
Multi-band accel. factor	3

**Geometry - AutoAlign**

Slice group	1
Position	R1.9 P25.7 H32.9 mm
Orientation	T > C-15.4
Phase enc. dir.	A >> P
AutoAlign	---
Initial Position	R1.9 P25.7 H32.9
R	1.9 mm
P	25.7 mm
H	32.9 mm
Initial Rotation	0.00 deg
Initial Orientation	T > C
T > C	-15.4
> S	0.0

**Geometry - Saturation**

Fat suppr.	Fat sat.
Special sat.	None

**System - Miscellaneous**

Positioning mode	FIX
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Sum of Squares
Matrix Optimization	Off
AutoAlign	---
Coil Select Mode	Default

**System - Adjustments**

B0 Shim mode	Standard
B1 Shim mode	TrueForm
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Dominant Fat	Off

**System - Adjustments**

Assume Silicone	Off
Adjustment Tolerance	Auto

**System - Adjust Volume**

Position	R1.9 P25.7 H32.9 mm
Orientation	T > C-15.4
Rotation	0.00 deg
A >> P	204 mm
R >> L	204 mm
F >> H	127 mm
Reset	Off

**System - Tx/Rx**

Frequency 1H	123.253584 MHz
Correction factor	1
Gain	High
Img. Scale Cor.	1.000
Reset	Off
? Ref. amplitude 1H	0.000 V

**Physio - Signal1**

1st Signal/Mode	None
TR	2000 ms
Multi-band accel. factor	3

**BOLD**

GLM Statistics	Off
Dynamic t-maps	Off
Ignore meas. at start	0
Ignore after transition	0
Model transition states	On
Temp. highpass filter	On
Threshold	4.00
Paradigm size	20
Meas[1]	Baseline
Meas[2]	Baseline
Meas[3]	Baseline
Meas[4]	Baseline
Meas[5]	Baseline
Meas[6]	Baseline
Meas[7]	Baseline
Meas[8]	Baseline
Meas[9]	Baseline
Meas[10]	Baseline
Meas[11]	Active
Meas[12]	Active
Meas[13]	Active
Meas[14]	Active
Meas[15]	Active
Meas[16]	Active
Meas[17]	Active
Meas[18]	Active
Meas[19]	Active
Meas[20]	Active
Motion correction	Off
Spatial filter	Off
Measurements	178
Delay in TR	0 ms
Multiple series	Off

**Sequence - Part 1**

Introduction	Off
Contrasts	1
Flow comp.	No

**Sequence - Part 1**

Multi-slice mode	Interleaved
Free echo spacing	Off
Echo spacing	0.77 ms
Bandwidth	1552 Hz/Px

**Sequence - Part 2**

EPI factor	92
Gradient mode	Normal
RF spoiling	Off

**Sequence - Special**

Excite pulse duration	4380 us
Single-band images	Off
MB LeakBlock kernel	On
MB dual kernel	Off
MB RF phase scramble	Off
SENSE1 coil combine	On
Invert RO/PE polarity	Off
PF omits higher k-space	Off
Disable freq. update	Off
Force equal slice timing	Off
Online multi-band recon.	Online
FFT scale factor	1.00
Physio recording	Off
Triggering scheme	Standard

\\USER\AP\TEST\EMPRISE\_220217\cmrr\_mbep2d\_bold\_2.2iso\_mod270122

TA: 6:06 PM: FIX Voxel size: 2.2×2.2×2.2 mmPAT: Off Rel. SNR: 1.00 : epfid

**Properties**

Prio recon	Off
Load images to viewer	On
Inline movie	Off
Auto store images	On
Load images to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Auto close inline display	Off
Start measurement without further preparation	Off
Wait for user to start	On
Start measurements	Single measurement

**Routine**

Slice group	1
Slices	51
Dist. factor	13 %
Position	R1.9 P25.7 H32.9 mm
Orientation	T > C-15.4
Phase enc. dir.	A >> P
AutoAlign	---
Phase oversampling	0 %
FoV read	204 mm
FoV phase	100.0 %
Slice thickness	2.20 mm
TR	2000 ms
TE	24.60 ms
Multi-band accel. factor	3
Filter	Prescan Normalize
Coil elements	HEA;HEP

**Contrast - Common**

TR	2000 ms
TE	24.60 ms
MTC	Off
Magn. preparation	None
Flip angle	80 deg
Fat suppr.	Fat sat.

**Contrast - Dynamic**

Averaging mode	Long term
Reconstruction	Magnitude
Measurements	178
Delay in TR	0 ms
Multiple series	Off

**Resolution - Common**

FoV read	204 mm
FoV phase	100.0 %
Slice thickness	2.20 mm
Base resolution	92
Phase resolution	100 %
Phase partial Fourier	6/8
Interpolation	Off

**Resolution - iPAT**

PAT mode	None
----------	------

**Resolution - Filter Image**

Distortion Corr.	Off
------------------	-----

**Resolution - Filter Image**

Prescan Normalize	On
-------------------	----

**Resolution - Filter Rawdata**

Raw filter	Off
Elliptical filter	Off
Hamming	Off

**Geometry - Common**

Slice group	1
Slices	51
Dist. factor	13 %
Position	R1.9 P25.7 H32.9 mm
Orientation	T > C-15.4
Phase enc. dir.	A >> P
FoV read	204 mm
FoV phase	100.0 %
Slice thickness	2.20 mm
TR	2000 ms
Multi-slice mode	Interleaved
Series	Interleaved
Multi-band accel. factor	3

**Geometry - AutoAlign**

Slice group	1
Position	R1.9 P25.7 H32.9 mm
Orientation	T > C-15.4
Phase enc. dir.	A >> P
AutoAlign	---
Initial Position	R1.9 P25.7 H32.9
R	1.9 mm
P	25.7 mm
H	32.9 mm
Initial Rotation	0.00 deg
Initial Orientation	T > C
T > C	-15.4
> S	0.0

**Geometry - Saturation**

Fat suppr.	Fat sat.
Special sat.	None

**System - Miscellaneous**

Positioning mode	FIX
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Sum of Squares
Matrix Optimization	Off
AutoAlign	---
Coil Select Mode	Default

**System - Adjustments**

B0 Shim mode	Standard
B1 Shim mode	TrueForm
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Dominant Fat	Off

**System - Adjustments**

Assume Silicone	Off
Adjustment Tolerance	Auto

**System - Adjust Volume**

Position	R1.9 P25.7 H32.9 mm
Orientation	T > C-15.4
Rotation	0.00 deg
A >> P	204 mm
R >> L	204 mm
F >> H	127 mm
Reset	Off

**System - Tx/Rx**

Frequency 1H	123.253584 MHz
Correction factor	1
Gain	High
Img. Scale Cor.	1.000
Reset	Off
? Ref. amplitude 1H	0.000 V

**Physio - Signal1**

1st Signal/Mode	None
TR	2000 ms
Multi-band accel. factor	3

**BOLD**

GLM Statistics	Off
Dynamic t-maps	Off
Ignore meas. at start	0
Ignore after transition	0
Model transition states	On
Temp. highpass filter	On
Threshold	4.00
Paradigm size	20
Meas[1]	Baseline
Meas[2]	Baseline
Meas[3]	Baseline
Meas[4]	Baseline
Meas[5]	Baseline
Meas[6]	Baseline
Meas[7]	Baseline
Meas[8]	Baseline
Meas[9]	Baseline
Meas[10]	Baseline
Meas[11]	Active
Meas[12]	Active
Meas[13]	Active
Meas[14]	Active
Meas[15]	Active
Meas[16]	Active
Meas[17]	Active
Meas[18]	Active
Meas[19]	Active
Meas[20]	Active
Motion correction	Off
Spatial filter	Off
Measurements	178
Delay in TR	0 ms
Multiple series	Off

**Sequence - Part 1**

Introduction	Off
Contrasts	1
Flow comp.	No

**Sequence - Part 1**

Multi-slice mode	Interleaved
Free echo spacing	Off
Echo spacing	0.77 ms
Bandwidth	1552 Hz/Px

**Sequence - Part 2**

EPI factor	92
Gradient mode	Normal
RF spoiling	Off

**Sequence - Special**

Excite pulse duration	4380 us
Single-band images	Off
MB LeakBlock kernel	On
MB dual kernel	Off
MB RF phase scramble	Off
SENSE1 coil combine	On
Invert RO/PE polarity	Off
PF omits higher k-space	Off
Disable freq. update	Off
Force equal slice timing	Off
Online multi-band recon.	Online
FFT scale factor	1.00
Physio recording	Off
Triggering scheme	Standard

\\USER\AP\TEST\EMPRISE\_220217\cmrr\_mbep2d\_bold\_2.2iso\_mod270122

TA: 6:06 PM: FIX Voxel size: 2.2×2.2×2.2 mmPAT: Off Rel. SNR: 1.00 : epfid

**Properties**

Prio recon	Off
Load images to viewer	On
Inline movie	Off
Auto store images	On
Load images to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Auto close inline display	Off
Start measurement without further preparation	Off
Wait for user to start	On
Start measurements	Single measurement

**Routine**

Slice group	1
Slices	51
Dist. factor	13 %
Position	R1.9 P25.7 H32.9 mm
Orientation	T > C-15.4
Phase enc. dir.	A >> P
AutoAlign	---
Phase oversampling	0 %
FoV read	204 mm
FoV phase	100.0 %
Slice thickness	2.20 mm
TR	2000 ms
TE	24.60 ms
Multi-band accel. factor	3
Filter	Prescan Normalize
Coil elements	HEA;HEP

**Contrast - Common**

TR	2000 ms
TE	24.60 ms
MTC	Off
Magn. preparation	None
Flip angle	80 deg
Fat suppr.	Fat sat.

**Contrast - Dynamic**

Averaging mode	Long term
Reconstruction	Magnitude
Measurements	178
Delay in TR	0 ms
Multiple series	Off

**Resolution - Common**

FoV read	204 mm
FoV phase	100.0 %
Slice thickness	2.20 mm
Base resolution	92
Phase resolution	100 %
Phase partial Fourier	6/8
Interpolation	Off

**Resolution - iPAT**

PAT mode	None
----------	------

**Resolution - Filter Image**

Distortion Corr.	Off
------------------	-----

**Resolution - Filter Image**

Prescan Normalize	On
-------------------	----

**Resolution - Filter Rawdata**

Raw filter	Off
Elliptical filter	Off
Hamming	Off

**Geometry - Common**

Slice group	1
Slices	51
Dist. factor	13 %
Position	R1.9 P25.7 H32.9 mm
Orientation	T > C-15.4
Phase enc. dir.	A >> P
FoV read	204 mm
FoV phase	100.0 %
Slice thickness	2.20 mm
TR	2000 ms
Multi-slice mode	Interleaved
Series	Interleaved
Multi-band accel. factor	3

**Geometry - AutoAlign**

Slice group	1
Position	R1.9 P25.7 H32.9 mm
Orientation	T > C-15.4
Phase enc. dir.	A >> P
AutoAlign	---
Initial Position	R1.9 P25.7 H32.9
R	1.9 mm
P	25.7 mm
H	32.9 mm
Initial Rotation	0.00 deg
Initial Orientation	T > C
T > C	-15.4
> S	0.0

**Geometry - Saturation**

Fat suppr.	Fat sat.
Special sat.	None

**System - Miscellaneous**

Positioning mode	FIX
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Sum of Squares
Matrix Optimization	Off
AutoAlign	---
Coil Select Mode	Default

**System - Adjustments**

B0 Shim mode	Standard
B1 Shim mode	TrueForm
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Dominant Fat	Off

**System - Adjustments**

Assume Silicone	Off
Adjustment Tolerance	Auto

**System - Adjust Volume**

Position	R1.9 P25.7 H32.9 mm
Orientation	T > C-15.4
Rotation	0.00 deg
A >> P	204 mm
R >> L	204 mm
F >> H	127 mm
Reset	Off

**System - Tx/Rx**

Frequency 1H	123.253584 MHz
Correction factor	1
Gain	High
Img. Scale Cor.	1.000
Reset	Off
? Ref. amplitude 1H	0.000 V

**Physio - Signal1**

1st Signal/Mode	None
TR	2000 ms
Multi-band accel. factor	3

**BOLD**

GLM Statistics	Off
Dynamic t-maps	Off
Ignore meas. at start	0
Ignore after transition	0
Model transition states	On
Temp. highpass filter	On
Threshold	4.00
Paradigm size	20
Meas[1]	Baseline
Meas[2]	Baseline
Meas[3]	Baseline
Meas[4]	Baseline
Meas[5]	Baseline
Meas[6]	Baseline
Meas[7]	Baseline
Meas[8]	Baseline
Meas[9]	Baseline
Meas[10]	Baseline
Meas[11]	Active
Meas[12]	Active
Meas[13]	Active
Meas[14]	Active
Meas[15]	Active
Meas[16]	Active
Meas[17]	Active
Meas[18]	Active
Meas[19]	Active
Meas[20]	Active
Motion correction	Off
Spatial filter	Off
Measurements	178
Delay in TR	0 ms
Multiple series	Off

**Sequence - Part 1**

Introduction	Off
Contrasts	1
Flow comp.	No

**Sequence - Part 1**

Multi-slice mode	Interleaved
Free echo spacing	Off
Echo spacing	0.77 ms
Bandwidth	1552 Hz/Px

**Sequence - Part 2**

EPI factor	92
Gradient mode	Normal
RF spoiling	Off

**Sequence - Special**

Excite pulse duration	4380 us
Single-band images	Off
MB LeakBlock kernel	On
MB dual kernel	Off
MB RF phase scramble	Off
SENSE1 coil combine	On
Invert RO/PE polarity	Off
PF omits higher k-space	Off
Disable freq. update	Off
Force equal slice timing	Off
Online multi-band recon.	Online
FFT scale factor	1.00
Physio recording	Off
Triggering scheme	Standard

\\USER\AP\TEST\EMPRISE\_220217\cmrr\_mbep2d\_bold\_2.2iso\_mod270122\_loc\_word

TA: 1:10 PM: FIX Voxel size: 2.2×2.2×2.2 mmPAT: Off Rel. SNR: 1.00 : efpid

**Properties**

Prio recon	Off
Load images to viewer	On
Inline movie	Off
Auto store images	On
Load images to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Auto close inline display	Off
Start measurement without further preparation	Off
Wait for user to start	On
Start measurements	Single measurement

**Routine**

Slice group	1
Slices	51
Dist. factor	13 %
Position	R1.9 P25.7 H32.9 mm
Orientation	T > C-15.4
Phase enc. dir.	A >> P
AutoAlign	---
Phase oversampling	0 %
FoV read	204 mm
FoV phase	100.0 %
Slice thickness	2.20 mm
TR	2000 ms
TE	24.60 ms
Multi-band accel. factor	3
Filter	Prescan Normalize
Coil elements	HEA;HEP

**Contrast - Common**

TR	2000 ms
TE	24.60 ms
MTC	Off
Magn. preparation	None
Flip angle	80 deg
Fat suppr.	Fat sat.

**Contrast - Dynamic**

Averaging mode	Long term
Reconstruction	Magnitude
Measurements	30
Delay in TR	0 ms
Multiple series	Off

**Resolution - Common**

FoV read	204 mm
FoV phase	100.0 %
Slice thickness	2.20 mm
Base resolution	92
Phase resolution	100 %
Phase partial Fourier	6/8
Interpolation	Off

**Resolution - iPAT**

PAT mode	None
----------	------

**Resolution - Filter Image**

Distortion Corr.	Off
------------------	-----

**Resolution - Filter Image**

Prescan Normalize	On
-------------------	----

**Resolution - Filter Rawdata**

Raw filter	Off
Elliptical filter	Off
Hamming	Off

**Geometry - Common**

Slice group	1
Slices	51
Dist. factor	13 %
Position	R1.9 P25.7 H32.9 mm
Orientation	T > C-15.4
Phase enc. dir.	A >> P
FoV read	204 mm
FoV phase	100.0 %
Slice thickness	2.20 mm
TR	2000 ms
Multi-slice mode	Interleaved
Series	Interleaved
Multi-band accel. factor	3

**Geometry - AutoAlign**

Slice group	1
Position	R1.9 P25.7 H32.9 mm
Orientation	T > C-15.4
Phase enc. dir.	A >> P
AutoAlign	---
Initial Position	R1.9 P25.7 H32.9
R	1.9 mm
P	25.7 mm
H	32.9 mm
Initial Rotation	0.00 deg
Initial Orientation	T > C
T > C	-15.4
> S	0.0

**Geometry - Saturation**

Fat suppr.	Fat sat.
Special sat.	None

**System - Miscellaneous**

Positioning mode	FIX
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Sum of Squares
Matrix Optimization	Off
AutoAlign	---
Coil Select Mode	Default

**System - Adjustments**

B0 Shim mode	Standard
B1 Shim mode	TrueForm
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Dominant Fat	Off



**System - Adjustments**

Assume Silicone	Off
Adjustment Tolerance	Auto

**System - Adjust Volume**

Position	R1.9 P25.7 H32.9 mm
Orientation	T > C-15.4
Rotation	0.00 deg
A >> P	204 mm
R >> L	204 mm
F >> H	127 mm
Reset	Off

**System - Tx/Rx**

Frequency 1H	123.253584 MHz
Correction factor	1
Gain	High
Img. Scale Cor.	1.000
Reset	Off
? Ref. amplitude 1H	0.000 V

**Physio - Signal1**

1st Signal/Mode	None
TR	2000 ms
Multi-band accel. factor	3

**BOLD**

GLM Statistics	Off
Dynamic t-maps	Off
Ignore meas. at start	0
Ignore after transition	0
Model transition states	On
Temp. highpass filter	On
Threshold	4.00
Paradigm size	20
Meas[1]	Baseline
Meas[2]	Baseline
Meas[3]	Baseline
Meas[4]	Baseline
Meas[5]	Baseline
Meas[6]	Baseline
Meas[7]	Baseline
Meas[8]	Baseline
Meas[9]	Baseline
Meas[10]	Baseline
Meas[11]	Active
Meas[12]	Active
Meas[13]	Active
Meas[14]	Active
Meas[15]	Active
Meas[16]	Active
Meas[17]	Active
Meas[18]	Active
Meas[19]	Active
Meas[20]	Active
Motion correction	Off
Spatial filter	Off
Measurements	30
Delay in TR	0 ms
Multiple series	Off

**Sequence - Part 1**

Introduction	Off
Contrasts	1
Flow comp.	No

**Sequence - Part 1**

Multi-slice mode	Interleaved
Free echo spacing	Off
Echo spacing	0.77 ms
Bandwidth	1552 Hz/Px

**Sequence - Part 2**

EPI factor	92
Gradient mode	Normal
RF spoiling	Off

**Sequence - Special**

Excite pulse duration	4380 us
Single-band images	Off
MB LeakBlock kernel	On
MB dual kernel	Off
MB RF phase scramble	Off
SENSE1 coil combine	On
Invert RO/PE polarity	Off
PF omits higher k-space	Off
Disable freq. update	Off
Force equal slice timing	Off
Online multi-band recon.	Online
FFT scale factor	1.00
Physio recording	Off
Triggering scheme	Standard

\\USER\AP\TEST\EMPRISE\_220217\cmrr\_mbep2d\_bold\_2.2iso\_mod270122\_loc\_visual

TA: 2:10 PM: FIX Voxel size: 2.2×2.2×2.2 mmPAT: Off Rel. SNR: 1.00 : epfid

**Properties**

Prio recon	Off
Load images to viewer	On
Inline movie	Off
Auto store images	On
Load images to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Auto close inline display	Off
Start measurement without further preparation	Off
Wait for user to start	On
Start measurements	Single measurement

**Routine**

Slice group	1
Slices	51
Dist. factor	13 %
Position	R1.9 P25.7 H32.9 mm
Orientation	T > C-15.4
Phase enc. dir.	A >> P
AutoAlign	---
Phase oversampling	0 %
FoV read	204 mm
FoV phase	100.0 %
Slice thickness	2.20 mm
TR	2000 ms
TE	24.60 ms
Multi-band accel. factor	3
Filter	Prescan Normalize
Coil elements	HEA;HEP

**Contrast - Common**

TR	2000 ms
TE	24.60 ms
MTC	Off
Magn. preparation	None
Flip angle	80 deg
Fat suppr.	Fat sat.

**Contrast - Dynamic**

Averaging mode	Long term
Reconstruction	Magnitude
Measurements	60
Delay in TR	0 ms
Multiple series	Off

**Resolution - Common**

FoV read	204 mm
FoV phase	100.0 %
Slice thickness	2.20 mm
Base resolution	92
Phase resolution	100 %
Phase partial Fourier	6/8
Interpolation	Off

**Resolution - iPAT**

PAT mode	None
----------	------

**Resolution - Filter Image**

Distortion Corr.	Off
------------------	-----

**Resolution - Filter Image**

Prescan Normalize	On
-------------------	----

**Resolution - Filter Rawdata**

Raw filter	Off
Elliptical filter	Off
Hamming	Off

**Geometry - Common**

Slice group	1
Slices	51
Dist. factor	13 %
Position	R1.9 P25.7 H32.9 mm
Orientation	T > C-15.4
Phase enc. dir.	A >> P
FoV read	204 mm
FoV phase	100.0 %
Slice thickness	2.20 mm
TR	2000 ms
Multi-slice mode	Interleaved
Series	Interleaved
Multi-band accel. factor	3

**Geometry - AutoAlign**

Slice group	1
Position	R1.9 P25.7 H32.9 mm
Orientation	T > C-15.4
Phase enc. dir.	A >> P
AutoAlign	---
Initial Position	R1.9 P25.7 H32.9
R	1.9 mm
P	25.7 mm
H	32.9 mm
Initial Rotation	0.00 deg
Initial Orientation	T > C
T > C	-15.4
> S	0.0

**Geometry - Saturation**

Fat suppr.	Fat sat.
Special sat.	None

**System - Miscellaneous**

Positioning mode	FIX
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Sum of Squares
Matrix Optimization	Off
AutoAlign	---
Coil Select Mode	Default

**System - Adjustments**

B0 Shim mode	Standard
B1 Shim mode	TrueForm
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Dominant Fat	Off

**System - Adjustments**

Assume Silicone	Off
Adjustment Tolerance	Auto

**System - Adjust Volume**

Position	R1.9 P25.7 H32.9 mm
Orientation	T > C-15.4
Rotation	0.00 deg
A >> P	204 mm
R >> L	204 mm
F >> H	127 mm
Reset	Off

**System - Tx/Rx**

Frequency 1H	123.253584 MHz
Correction factor	1
Gain	High
Img. Scale Cor.	1.000
Reset	Off
? Ref. amplitude 1H	0.000 V

**Physio - Signal1**

1st Signal/Mode	None
TR	2000 ms
Multi-band accel. factor	3

**BOLD**

GLM Statistics	Off
Dynamic t-maps	Off
Ignore meas. at start	0
Ignore after transition	0
Model transition states	On
Temp. highpass filter	On
Threshold	4.00
Paradigm size	20
Meas[1]	Baseline
Meas[2]	Baseline
Meas[3]	Baseline
Meas[4]	Baseline
Meas[5]	Baseline
Meas[6]	Baseline
Meas[7]	Baseline
Meas[8]	Baseline
Meas[9]	Baseline
Meas[10]	Baseline
Meas[11]	Active
Meas[12]	Active
Meas[13]	Active
Meas[14]	Active
Meas[15]	Active
Meas[16]	Active
Meas[17]	Active
Meas[18]	Active
Meas[19]	Active
Meas[20]	Active
Motion correction	Off
Spatial filter	Off
Measurements	60
Delay in TR	0 ms
Multiple series	Off

**Sequence - Part 1**

Introduction	Off
Contrasts	1
Flow comp.	No

**Sequence - Part 1**

Multi-slice mode	Interleaved
Free echo spacing	Off
Echo spacing	0.77 ms
Bandwidth	1552 Hz/Px

**Sequence - Part 2**

EPI factor	92
Gradient mode	Normal
RF spoiling	Off

**Sequence - Special**

Excite pulse duration	4380 us
Single-band images	Off
MB LeakBlock kernel	On
MB dual kernel	Off
MB RF phase scramble	Off
SENSE1 coil combine	On
Invert RO/PE polarity	Off
PF omits higher k-space	Off
Disable freq. update	Off
Force equal slice timing	Off
Online multi-band recon.	Online
FFT scale factor	1.00
Physio recording	Off
Triggering scheme	Standard